COVID N95 DECON & REUSE



DECON OR DISPOSE?





No soil
Nosepiece intact
Straps intact
No tears
No deformation







Soiling

Deformation



Broken Strap

Tearing





Nosepiece Detachment

Nosepiece Fracture



The Content provided by N95DECON is for INFORMATIONAL PURPOSES ONLY and DOES NOT CONSTITUTE THE PROVIDING OF MEDICAL ADVICE and IS NOT INTENDED TO BE A SUBSTITUTE FOR INDEPENDENT PROFESSIONAL MEDICAL JUDGMENT, ADVICE, DIAGNOSIS, OR TREATMENT. Use or reliance on any Content provided by N95DECON is SOLELY AT YOUR OWN RISK. A link to the full N95DECON disclaimer can be found at https://www.n95decon.org/disclaimer.



v1.0 (April 17, 2020) www.n95decon.org