



循環動態維持を補助する システムをキット化。 セットアップもスピーディです。

経皮的循環補助システム

キャピオックスEBS



Quick

- 遠心ポンプ、人工肺、血液回路をキット化。短時間で簡単にセットアップ可能。
- オートプライミング機能により、約5分でプライミング操作を完了。
- ロックコネクター式のカテーテルで、回路との脱着が簡単。

Compact

- プライミングボリュームは470mlと低減化を実現。

Simple

- 熱交換器やエア抜き用リザーバー、フィルター等を省いたシンプルな回路構成。

Safety

- 停電時は内蔵バッテリーに自動的に切り換え。
- 血液損傷の少ない直流流路構造の遠心ポンプ。



承認番号：キャピオックスEBS 心肺キット [6B-207] キャピオックスEBS カテーテルキット [4B-421]
キャピオックス遠心ポンプコントローラーSP101 [6B-159]

テルモ株式会社 〒151 東京都渋谷区幡ヶ谷2-44-1 ④、TERUMO、キャピオックス、CAPIOXは、テルモ株式会社の登録商標です。

エマセブ、EMERSAVEは、テルモ株式会社の商標です。




CAPIOX[®] SX 10 designed for small-volume requirements.

- 135ml* low priming volume
- Extra-small dynamic priming volume
- Efficient blood volume management
- Outstanding gas exchange performance
- Superior heat exchange system
- Open or closed system flexibility

* For larger volume requirements, Terumo's standard 270ml CAPIOX SX 18 oxygenator system is available.

135ml
Priming
Volume



CAPIOX[®] SX 10

Hollow Fiber Oxygenator/Detachable Hardshell Reservoir
for Pediatric and Small-Volume Requirements

 **TERUMO[®]**

TERUMO CORPORATION 44-1, 2-chome, Hatagaya, Shibuya-ku, Tokyo 151, Japan Tel: 81-3-3374-8111 Fax: 81-3-3374-8196
キャビオックス SX [5B-61]  and TERUMO are registered trademarks of TERUMO CORPORATION.